


SOL-M™ Blunt Fill Needles & Blunt Fill Needle with 5 micron filter

SOL-M™ Blunt Fill Needles are used for draw-up of medication, re-constitution of medication and accessing IV Bags.

SOL-M™ Blunt Fill Needle with filter features a 5 micron filter membrane that is designed to prevent larger particles being drawn into the syringe when preparing medications from ampules.

Features & Benefits

- **45°** bevel for ease in draw-up of medication from vial
 - **18G** for ease in draw-up of medications
 - **Conventional standard hub** to fit standard luer lock syringes
 - **5 micron** filter - Effective filtration of contaminant particles.
- Color coded for identification
Lower cost than sharps engineered devices



REF	Description	Qty
110022	■ 18G x 1 1/2"	100/Box; 1000/Case
110022F	■ 18G x 1 1/2" with 5 micron filter	100/Box; 1000/Case
110023	■ 18G x 50mm	100/Box; 1000/Case
110023F	■ 18G x 50mm with 5 micron filter	100/Box; 1000/Case

SOL-MILLENNIUM®

Shanghai Office:
Suite 701, Block E, 18 Dongfang Rd
Pudong, Shanghai, China, 200120
Tel: (86 21) 6862 2586
Fax: (86 21) 6862 2585
E mail: info@sol-millennium.com

North America Office:
5415 Sugarloaf Parkway, Suite 2203
Lawrenceville, GA 30043
Tel: 404.973.2200
Fax: 404.973.2212
E mail: info@sol-millenniumna.com

Europe Office:
Emilia Plater 53, PL-00113 Warszawa
Tel: +48 22 528 69 41
Fax: +48 22 528 67 01
Email: info@sol-millenniumeu.com

Australia Office:
Suite 404/8-12 Station St,
Homebush, NSW, 2140
Tel: 02 9746 7738
E-mail: vera@sol-millennium.com.au

Brasil Office:
Avenida Anhanguera, 869
Ribeirão Preto - São Paulo, Brasil
Tel: 14.025-480
Telefone: 55 16 3442 1454 / 55 11 984060126
E-mail: RBartolo@sol-millenniumbr.com