

TWIST UP&DOWN REFILL



Feature

- Easily removable metal spring
 - By one hand from elevating, pumping to closing.
 - Turn the shoulder and press the actuator to use it.
 - Up & down bottle is made PET for easy disposal.
-

How to refill



Ready inner bottle
& outer bottle

Put inner bottle
into outer bottle

① Turn Shoulder
② Inner bottle down

Finish

Available Options

Item



Item No.

FS TUD M / FS TUD E

Capacity

30/50 ml

Size

Ø 39 X 111.7 mm

Dosage

0.2cc / 0.08cc

Material

Inner Bottle : PP / PCR PP
Actuator & Nozzle : PP / PCR PP
Pump Cap : PP / PCR PP

Outer Bottle : PET / PCR PET
Shoulder : PET / PCR PET
Lifter : PET / PCR PET

Refill cap : PP
Stopper ring : PP

Button Type



Mist Pump



Essence Pump