


ARIANE

PUMP PACK : 2 FORMULATIONS TO MIX



- PP / PE / ELASTOMER

Ø20 - Ariane Pump

15 ml
H: 157 mm




ZOOM ON THE PRODUCT

A system that give the feeling to the consumer to take part in the ultimate process

Available in cream or spray version, dose from 130 µl or 230 µl

Packaging delivered in 3 parts (a bole with a pump, a tank and a valve)

Product distribution with a viscosity from 1 to 30000 centipoises



The Ariane Pack allows two formulations to be activated at the time of use:

- 1 REMOVE tamper tab
- 2 HIT the tank to release the active product: powder or liquid
- 3 SHACK the pack
- 4 APPLY pressure on the pump to distribute the product



1



2



3



4

DECORATION

Screen (silk) printing

Offset

Hot stamping

You can create your own artwork without technical limit.

Ø20 ARIANE PUMP

DISTRIBUTION : Liquids, Powders, Serum, Cream